

Product	:	Anti-Testosterone
Cat. No.	:	PAS9247
Lot No.	:	RS568-08
Immunogen	:	17 β Testosterone (7 α)-BTG
Host Animal	:	Sheep
Format	:	Ig fraction
Ig concentration	:	6.94mg/ml (U.V.abs @ 280nm)
Working Concentration	:	1.25 μ g/ml (EIA)
PH	:	7.2
Buffer	:	20mM Phosphate, 150mM Sodium Chloride
Preservative	:	0.09% Sodium Azide
Storage	:	-20°C