

Product	:	Anti-Leucomalachite Green
Cat. No.	:	PAS9533
Lot No.	:	RS476-07/BC4
Immunogen	:	Leucomalachite Green-BTG
Host Animal	:	Sheep
Format	:	Ig fraction
Ig Concentration	:	8.58mg/ml (U.V. abs @ 280nm)
Working Concentration	:	0.625µg/ml (EIA)
pH	:	7.2
Buffer	:	20mM Phosphate, 150mM Sodium Chloride
Preservative	:	0.09% Sodium Azide
Storage	:	-20°C

Specificity

Compound	% Cross-Reactivity
Leucomalachite Green	100
Malachite Green	26.4
Leucocrystal Violet	200