

<b>Product</b>	:	Anti-17 $\alpha$ Nortestosterone
<b>Cat. No.</b>	:	PAS9556
<b>Lot No.</b>	:	RS599-08
<b>Immunogen</b>	:	17 $\alpha$ -19- Nortestosterone(3)-BTG
<b>Host Animal</b>	:	Sheep
<b>Format</b>	:	Ig fraction
<b>Ig Concentration</b>	:	7.56mg/ml (U.V.abs @ 280nm)
<b>Working Concentration</b>	:	0.25 $\mu$ g/ml (EIA)
<b>PH</b>	:	7.2
<b>Buffer</b>	:	20mM Phosphate, 150mM Sodium Chloride
<b>Preservative</b>	:	0.09% Sodium Azide
<b>Storage</b>	:	-20 $^{\circ}$ C