

Product	:	Anti-17 α - Boldenone
Cat. No.	:	PAS9677
Lot No.	:	RS686-10
Immunogen	:	17 α - Boldenone (3)-BTG
Host Animal	:	Sheep
Format	:	Ig fraction
Ig Concentration	:	13.80mg/ml (U.V. abs @ 280nm)
Working Concentration	:	0.25 μ g/ml (EIA)
PH	:	7.2
Buffer	:	20mM Phosphate, 150mM Sodium Chloride
Preservative	:	0.09% Sodium Azide
Storage	:	-20°C