

Product	:	Anti-Creatinine
Cat. No.	:	PAS 9686
Lot No.	:	RS300-24
Immunogen	:	Creatinine (N)-BSA
Host Animal	:	Sheep
Format	:	Ig fraction
Ig Concentration	:	6.13mg/ml (U.V. abs @ 280nm)
Working Concentration	:	2µg/ml (EIA)
PH	:	7.2
Buffer	:	20mM Phosphate, 150mM Sodium Chloride
Preservative	:	0.09% Sodium Azide
Storage	:	-20°C