

Product	:	Anti-Fentanyl/Despropionyl Fentanyl
Cat. No.	:	PAS 9694
Lot No.	:	RS167-12
Immunogen	:	Fentanyl (p)-BSA
Host Animal	:	Sheep
Format	:	Ig Fraction
Ig Concentration	:	7.67mg/ml (U.V. abs @ 280nm)
Working Concentration	:	1.0µg/ml (EIA)
PH	:	7.2
Buffer	:	20mM Phosphate, 150mM Sodium Chloride
Preservative	:	0.09% Sodium Azide
Storage	:	-20°C