

Product	:	Anti-B-Lactam Penicillin
Cat. No.	:	PAS9718
Lot No.	:	RS16-25
Immunogen	:	Ampicillin-BSA
Host Animal	:	Sheep
Format	:	Ig Fraction
Ig Concentration	:	8.14mg/ml (U.V abs @ 280nm)
Working Concentration	:	8.14µg/ml (EIA)
PH	:	7.2
Buffer	:	20mM Phosphate, 150mM Sodium Chloride
Preservative	:	0.09% Sodium Azide
Storage	:	-20°C

Specificity

Compound	% Cross-Reactivity
Ampicillin	100
Amoxycillin	87
Penicillin G	72
Oxacillin	39
Cloxacillin	45
Dicloxacillin	46