

Product	:	Anti-Ractopamine
Cat. No.	:	PAS 9725
Lot No.	:	RS363-B15
Immunogen	:	Ractopamine (10)-BSA
Host Animal	:	Sheep
Format	:	Ig fraction
Ig Concentration	:	6.88mg/ml (U.V. abs @ 280nm)
Working Concentration	:	1.25µg/ml (EIA)
PH	:	7.2
Buffer	:	20mM Phosphate, 150mM Sodium Chloride
Preservative	:	0.09% Sodium Azide
Storage	:	-20°C

Specificity

Compound	% Cross-Reactivity
Ractopamine	100
Dobutamine	<0.8
Isoxsuprine	<0.6
DL-Isoproterenol	<0.5
Ritodrine	<0.5
Salmeterol Xinafoate	<0.5
Methyl Clenbuterol	<0.05
Fenoterol	<0.05
Pirbuterol	<0.05
Tulobuterol	<0.05
Dopamine	<0.01
Salbutamol	<0.01
Clenbuterol	<0.01
Metaproterenol	<0.01
Terbutaline	<0.01