

Product	:	Anti-Testosterone
Cat. No.	:	PAS9745
Lot No.	:	RS323-19
Immunogen	:	Testosterone (3) - long spacer - BSA
Host Animal	:	Sheep
Format	:	Ig fraction
Ig concentration	:	6.02mg/ml (U.V.abs @ 280nm)
Working Concentration	:	2.5µg/ml (EIA)
PH	:	7.2
Buffer	:	20mM Phosphate, 150mM Sodium Chloride
Preservative	:	0.09% Sodium Azide
Storage	:	-20°C