

Product	:	Anti-Mycophenolic Acid
Cat. No.	:	PAS9877
Lot No.	:	RS1260
Immunogen	:	Mycophenolic-BTG
Host Animal	:	Sheep
Format	:	Ig fraction
Ig concentration	:	8.80mg/ml (U.V.abs @ 280nm)
Working Concentration	:	0.625 µg/ml (EIA)
PH	:	7.2
Buffer	:	20mM Phosphate, 150mM Sodium Chloride
Preservative	:	0.09% Sodium Azide
Storage	:	-20°C