

PRODUCT SPECIFICATION

Product	:	Anti-Leucomalachite Green
Cat. No.	:	PAS 9805
Lot No.	:	RS476-07/BC8
Immunogen	:	Leucomalachite Green-BTG
Host Animal	:	Sheep
Format	:	Ig fraction
Ig Concentration	:	7.99mg/ml (U.V.abs @ 280nm)
Working Concentration	:	0.625µg/ml
pH	:	7.2
Buffer	:	20mM Phosphate, 150mM Sodium Chloride
Preservative	:	0.09% Sodium Azide
Storage	:	-20°C

SPECIFICITY:

Compound	% Cross-Reactivity
Leucomalachite Green	100
Malachite Green	26.4