

| | | |
|------------------------------|---|---------------------------------------|
| Product | : | Anti-Nitrofurantoin |
| Cat. No. | : | PAS9228 |
| Lot No. | : | RS494-12/BC1 |
| Immunogen | : | Aminohydantoin-BTG |
| Host Animal | : | Sheep |
| Format | : | Ig fraction |
| Ig concentration | : | 6.73mg/ml (U.V.abs @ 280nm) |
| Working Concentration | : | 5.0µg/ml (EIA) |
| PH | : | 7.2 |
| Buffer | : | 20mM Phosphate, 150mM Sodium Chloride |
| Preservative | : | 0.09% Sodium Azide |
| Storage | : | -20°C |