

Product	:	Anti-17 α Testosterone
Cat. No.	:	PAS9670
Lot No.	:	RS680-08
Immunogen	:	17 α Testosterone (3)-BTG
Host Animal	:	Sheep
Format	:	Ig fraction
Ig Concentration	:	9.84mg/ml (U.V. abs @ 280nm)
Working Concentration	:	2.5 μ g/ml (EIA)
pH	:	7.2
Buffer	:	20mM Phosphate, 150mM Sodium Chloride
Preservative	:	0.09% Sodium Azide
Storage	:	-20°C