

Product	:	Anti-17 α Nortestosterone
Cat. No.	:	PAS9707
Lot No.	:	RS599-08
Immunogen	:	17 α -19- Nortestosterone(3)-BTG
Host Animal	:	Sheep
Format	:	Ig fraction
Ig Concentration	:	7.56mg/ml (U.V.abs @ 280nm)
Working Concentration	:	0.25 μ g/ml (EIA)
PH	:	7.2
Buffer	:	20mM Phosphate, 150mM Sodium Chloride
Preservative	:	0.09% Sodium Azide
Storage	:	-20 $^{\circ}$ C