

Product	:	Anti-Oxazepam
Cat. No.	:	PAS9710
Lot No.	:	RS278-B12
Immunogen	:	Oxazepam (N)-BSA
Host Animal	:	Sheep
Format	:	Ig Fraction
Ig Concentration	:	5.84mg/ml (U.V. abs @ 280nm)
Working Concentration	:	1.25µg/ml (EIA)
PH	:	7.2
Buffer	:	20mM Phosphate, 150mM Sodium Chloride
Preservative	:	0.09% Sodium Azide