

Product	:	Anti-Oxazepam
Cat. No.	:	PAS9711
Lot No.	:	RS81-12
Immunogen	:	Oxazepam(3)-BSA
Host Animal	:	Sheep
Format	:	Ig Fraction
Ig Concentration	:	10.78mg/ml (U.V. abs @ 280nm)
Working Concentration	:	5.5µg/ml (EIA)
PH	:	7.2
Buffer	:	20mM Phosphate, 150mM Sodium Chloride
Preservative	:	0.09% Sodium Azide
Storage	:	-20°C

Specificity

Compound	% Cross-Reactivity
Oxazepam	100
Lorazepam	<0.77
Temazepam	<0.77
Alprazolam	<0.77
α-Hydroxyalprazolam	<0.77
Nordiazepam	974.68
Nitrazepam	55
7-Aminonitrazepam	<0.77
Clonazepam	<0.77