

<b>Product</b>	:	Anti-Testosterone
<b>Cat. No.</b>	:	PAS9745
<b>Lot No.</b>	:	RS323-19
<b>Immunogen</b>	:	Testosterone (3) - long spacer - BSA
<b>Host Animal</b>	:	Sheep
<b>Format</b>	:	Ig fraction
<b>Ig concentration</b>	:	6.02mg/ml (U.V.abs @ 280nm)
<b>Working Concentration</b>	:	2.5µg/ml (EIA)
<b>PH</b>	:	7.2
<b>Buffer</b>	:	20mM Phosphate, 150mM Sodium Chloride
<b>Preservative</b>	:	0.09% Sodium Azide
<b>Storage</b>	:	-20°C