

## CERTIFICATE OF ANALYSIS

<b>Product</b>	:	Ampicillin – HRP Conjugate
<b>Cat. No.</b>	:	HRP9339
<b>Lot No.</b>	:	ESC4149
<b>Volume</b>	:	1ml when reconstituted
<b>Preparation of Conjugate</b>	:	Add 1ml of Conjugate Diluent to 24ml of double deionised water. Reconstitute the contents of the vial of conjugate concentrate with 1ml of diluted conjugate diluent. Leave to stand for 10 minutes and then mix on a rocking roller for 10 minutes.
<b>Recommended Working Dilution (ELISA)</b>	:	1/100 when used in a competitive ELISA with antibody PAS9718. Dilute the reconstituted conjugate concentrate with diluted conjugate diluent.
<b>Storage</b>	:	Store the freeze-dried conjugate concentrate at 2 - 8°C. Conjugate concentrate is stable for 2 hours after reconstitution. Once the conjugate has been diluted, use immediately.