

Product:	Purified high affinity monoclonal antibody to salbutamol		
Immunogen:	Salbutamol-BTG		
Pack Sizes: Cat. No.	100 µg: MAB9345	Host Animal : Sheep	Lot No. 6B7.1508 08
	1 mg: MAB9901	Format : Purified monoclonal antibody	
Method of Purification: Protein A affinity chromatography			
Method of Production: In vitro cell culture using standard DMEM media supplemented with foetal bovine serum of a USDA approved origin.		Buffer: PBS at pH 7.3-7.4	
Concentration: 2.03 mg/ml (U.V. abs at 280nm with 1mg/ml = 1.4AU)		Preservative: 0.1% Sodium Azide	
Storage: +2-+8 ⁰ C.			
Purity Main peak ≥90% (by SEC-HPLC)			
Immunological activity Binds to salbutamol-BTG immunogen (ELISA method)			
Related Products: Salbutamol Polyclonal Antibodies: Cat. No PAS9727 size 100µl PAS9576 size 1ml Cat. No PAS9728 size 100µl PAS9577 size 1ml Tracer: Salbutamol-HRP Cat. No. HRP9348, size 0.5 ml			
This product as supplied is intended for research applications only, not for use in therapeutic or diagnostic applications without the expressed written authorization of Randox Life Sciences.			
Rev 1.0 April 10 Information correct at time of going to print			